

Galaxy Ultra 30D Mobile

Operating Theatre Lighting



- **Extended Arm Reach**
- **Ultra High Light Quality**
- **Excellent Manoeuvrability**
- **Emergency Battery Back-up**



Galaxy Ultra Mobile

Advanced Mobile Surgical Lighting

The latest extension to Brandon Medical's award winning Galaxy Ultra surgical light range. Offering the same ultra high performance and reliability as ceiling mounted lights but with infinite mobility and flexibility.



Extended Arm Reach

The unique design of the Galaxy Ultra Mobile base allows the light to be manoeuvred much closer to the table without impeding the surgeon. Small and compact, the reduced footprint still offers complete stability.

Ultra High Light Quality

A large diameter of illumination combines with ultra-high light intensity to provide cool, daylight quality light - giving accurate true-colour rendition in surgery. A unique reflector system allows "shadowless" performance, whilst focus control allows the beam to be adjusted without moving the light head.

Excellent Manoeuvrability

Easy to move on low inertia, anti-static, lockable double castor wheels. The light head is mounted on a spring-balanced arm for easy positioning.

Emergency Battery Back-up (option)

Located in the base, a sophisticated battery mechanism provides up to 3 hours of uninterrupted run-time in the event of power loss. The batteries are equipped with a "quick-charge" facility allowing ultra-fast battery charging.

Other features

- Ergonomic handle provides easy movement
- Cable storage keeps operating theatre tidy and free from trip hazards.
- Sterilisable handle allows surgeon to safely position light.
- Simple and user-friendly interface avoids any operational confusion.
- CE, BS EN 60601-1, ISO 9001:2000.



Brandon Medical – Brilliant by Design

Brandon Medical are experts in the field of medical lighting design with over 40 years of healthcare experience. The Company manufactures and markets a wide range of medical lighting products, from major operating theatre lights through to bed head lights. Brandon has a strong reputation for research and innovation, having won numerous international awards for design. Brandon products are in use throughout the world and are supported by an extensive distributor network.

It is the company's aim to deliver the best possible customer service, excellent quality and the most advanced products. Brandon Medical's policy is simply to be "Brilliant by Design".



**For more information
about Brandon products
call 0113 277 7393**

Technical data

Part number	Specifications
GU30DM	Galaxy Ultra 30D Mobile Operating Theatre Light
GU30DMEL	Galaxy Ultra 30D Mobile Operating Theatre Light with Emergency Battery Back-up

- **Light output** – 90,000 Lux @ 1m
- **Weight** – 95Kg (GU30DM)
– 102Kg (GU30DMEL)
- **Height** – 1620mm
- **Colour temperature** – 4,300°K
- **Colour Rendition Index R_a** – 96
- **Colour rendition R_g** – 92
- **Depth of light column** – 850mm
- **Infra-red absorption** – 99%
- **Focus range** – 700 – 1,600mm
from the lamp head
- **Battery life** – 3 hours
- **Lamp life** – 1,500 hours
- **Bulb** – 3 x 50W
- **Voltage S.E.L.V.** – 23 volts
- **Radiant energy units** – 3.86 W/m²/KLux
- **Handle sterilisation** – 134°C
- **Guarantee** – 3 years



Brandon Medical, Holme Well Road, Leeds LS10 4TQ
 Tel: +44 (0)113 2777 393 Fax: +44 (0)113 2728 844 Email: enquiries@brandon-medical.com
www.brandon-medical.com



In line with the Brandon policy of progressive improvement, the right is reserved to alter specifications and/or prices without notice